
Xanespy Documentation

Release 0.1

Mark Wolf

Jan 08, 2017

CONTENTS:

1	Introduction	3
1.1	X-Ray Absorbance Basics	3
1.2	Example Workflow	3
2	Importing Data into Xanespy	5
2.1	TXM from 8-BM (APS) and 6-2c (SSRL)	5
2.2	Ptychography from 5.3.2.1 (ALS)	5
3	Analyzing the Data	7
3.1	Frame Alignment	7
3.2	Subtracting Surroundings	7
3.3	Spectrum Fitting - K-Edge	7
3.4	Spectrum Fitting - L-Edge	7
4	Visualization of Results	9
4.1	Plotting	9
4.2	Interactive Viewer (Qt)	9
5	Accessing the Data Directly	11
6	xanespy	13
6.1	xanespy package	13
7	Indices and tables	45
	Python Module Index	47

Python library for analyzing **X-Ray absorbance spectroscopy** data.

Warning: This documentation is still under development. The code it represents could change at any time.

INTRODUCTION

Xanespy is a toolkit for interacting with X-ray microscopy data, most likely collected at a synchrotron beamline. By collecting a set of frames at multiple X-ray energies, spectral maps are reconstructed to provide chemical insight. Multiple framesets can be collected sequentially as part of an *operando* experiment and analyzed simultaneously in python. Slow operations take advantage of multiple cores when available.

1.1 X-Ray Absorbance Basics

1.2 Example Workflow

A typical procedure for interacting with microscope frame-sets involves the following parts:

- Import the raw data
- Apply corrections and align the images
- Calculate some metric and create maps of it
- Visualize the maps, static or interactively.

Example for a single frameset across an X-ray absorbance edge:

```
import xanespy

# Example for importing from SSRL beamline 6-2c
xanespy.import_ssrl_frameset('<data_dir>', hdf_filename='imported_data.h5')

# Load a pre-defined XAS edge or create your own subclass xanespy.Edge
edge = xanespy.k_edges['Ni_NCA']
# Now load the newly created HDF5 file and the X-ray absorbance edge
fs = xanespy.XanesFrameset(filename='imported_data.h5', edge=edge)

# Perform automatic frame alignment
fs.align_frames(passes=5)
# Fit the absorbance spectra and extract the edge position (SLOW!)
fs.fit_spectra()

# Inspect the result with the built-in Qt5 GUI
fs.qt_viewer()
```


IMPORTING DATA INTO XANESPY

2.1 TXM from 8-BM (APS) and 6-2c (SSRL)

2.2 Ptychography from 5.3.2.1 (ALS)

ANALYZING THE DATA

3.1 Frame Alignment

3.2 Subtracting Surroundings

3.3 Spectrum Fitting - K-Edge

3.4 Spectrum Fitting - L-Edge

VISUALIZATION OF RESULTS

4.1 Plotting

4.2 Interactive Viewer (Qt)

ACCESSING THE DATA DIRECTLY

While the `XanesFrameset` class has methods for common tasks, sometimes it is necessary to access the data directly, as either numpy arrays or h5py datasets. The `xanes_frameset` has a `store()` method that returns an interface (`TXMStore`) to the underlying HDF5 file.

Warning: The `TXMStore` created by `xanes_frameset.store()` is attached to an open HDF5 file. It is strongly recommended to use the `with` statement described below. Otherwise make sure to call the store's `close()` method in a `try...except` block. File corruption is possible if not opened in this manner.

Call the following to get access to the associated datasets. Properties of the interface will return an HDF5 dataset in most cases.:

```
import xanespy as xp
frameset = xp.XanesFrameset(...)

# Open the TXMStore interface
with frameset.store() as store:
    # For example, the images are in (timestep, energy, row, column) order
    assert store.absorbances.shape == (10, 62, 1024, 1024)
    # Energies are in (timestep, energy) order
    assert store.energies.shape == (10, 62)
```


XANESPY

6.1 xanespy package

6.1.1 Submodules

6.1.2 xanespy.beamlines module

Functions and classes that prepare experiments at specific synchrotron beamlines.

```
class xanespy.beamlines.Detector (start: xanespy.beamlines.ZoneplatePoint, z_step: int = None,  
                                end: xanespy.beamlines.ZoneplatePoint = None)
```

Bases: `xanespy.beamlines.Zoneplate`

A calibration object for the position of the detector.

```
class xanespy.beamlines.DetectorPoint (z, energy)
```

Bases: `tuple`

energy

Alias for field number 1

z

Alias for field number 0

```
class xanespy.beamlines.Zoneplate (start: xanespy.beamlines.ZoneplatePoint, z_step: int = None,  
                                end: xanespy.beamlines.ZoneplatePoint = None)
```

Bases: `object`

Type of focusing optic using in X-ray microscopy. It must be moved with changing energy to properly focus the beam. In order to properly predict zoneplate positions, it needs either two position-energy pairs or one position-energy pair and a step. Passing two position-energy pairs is preferred because this allows x, y and z to be set properly instead of just z.

Parameters

- **start** (–) –
- **z_step** (–) – of beam energy.
- **end** (–) –

position (*energy: float*)

Predict the x, y and z position of the zonplate for the given energy.

```
class xanespy.beamlines.ZoneplatePoint (x, y, z, energy)
```

Bases: `tuple`

energy

Alias for field number 3

x

Alias for field number 0

y

Alias for field number 1

z

Alias for field number 2

```
xanespy.beamlines.sector8_xanes_script (dest, edge: edges.KEdge, zoneplate: xane-  
spy.beamlines.Zoneplate, detector: xane-  
spy.beamlines.Detector, sample_positions: typ-  
ing.List[utilities.position], names: typing.List[str],  
iterations: typing.Iterable = range(0, 1), binning:  
int = 1, exposure: int = 30, abba_mode: bool =  
True)
```

Prepare an script file for running multiple consecutive XANES framesets on the transmission x-ray microscope at the Advanced Photon Source beamline 8-BM-B.

Parameters

- **dest** (-) –
- **edge** (-) –
- **binning** (-) – (eg. 2 means 2x2 CCD pixels become 1 image pixel.
- **exposure** (-) –
- **sample_positions** (-) – order to capture the image.
- **zoneplate** (-) –
- **detector** (-) –
- **names** (-) –
- **iterations** (-) – set of xanes location with reference.
- **abba_mode** (-) – to save time. Eg: reference, sample, change-energy, sample, reference, change-energy, etc. Not compatible with *frame_rest* argument.

```
xanespy.beamlines.ssrl6_xanes_script (dest, edge: edges.KEdge, zoneplate: xane-  
spy.beamlines.Zoneplate, positions: typ-  
ing.List[utilities.position], reference_position: util-  
ities.position, iterations: typing.Iterable, iteration_rest:  
int = 0, frame_rest: int = 0, binning: int = 2, exposure:  
int = 0.5, repetitions: int = 5, ref_repetitions: int = 10,  
abba_mode: bool = True)
```

Prepare a script file for running multiple consecutive XANES framesets on the transmission x-ray microscope at the Advanced Photon Source beamline 8-BM-B. Both *iteration_rest* and *frame_rest* can be used to give the material time to recover from X-ray damage.

Parameters

- **dest** (-) –
- **edge** (-) –
- **binning** (-) – (eg. 2 means 2x2 CCD pixels become 1 image pixel.
- **exposure** (-) –

- **positions** (-) – order to capture the image.
- **reference_position** (-) – reference frame.
- **iteration_rest** (-) – iterations. Beam will wait at reference location before starting next XANES set.
- **frame_rest** (-) – wait at reference location before starting next energy frame.
- **zoneplate** (-) –
- **detector** (-) –
- **iterations** (-) –
- **repetitions** (-) – location/energy. These frames will then be averaged during analysis.
- **ref_repetitions** (-) –
- **abba_mode** (-) – locations first to save time. Eg: reference, sample, change-energy, sample, reference, change-energy, etc. Not compatible with *frame_rest* argument.

`xanespy.beamlines.write_scaninfo_header` (*f*, *abba_mode*, *repetitions*, *ref_repetitions*)

6.1.3 xanespy.edges module

Descriptions of X-ray energy absorption edge.

class `xanespy.edges.Edge`

Bases: `object`

An X-ray absorption edge. It is defined by a series of energy ranges. All energies are assumed to be in units of electron-volts. This class is intended to be extended into K-edge, L-edge, etc.

E_0

float – The energy of the absorption edge itself.

regions

list of 3-tuples – All the energy regions. Each tuple is of the form (start, end, step) and is inclusive at both ends.

name

string – A human-readable name for this edge (eg “Ni K-edge”)

pre_edge

2-tuple – Energy range (start, stop) that defines points below the edge region, inclusive.

post_edge

2-tuple – Energy range (start, stop) that defines points above the edge region, inclusive.

post_edge_order

int – What degree polynomial to use for fitting the post_edge region.

map_range

2-tuple – Energy range (start, stop) used for normalizing maps. If not supplied, will be determine from pre- and post-edge arguments.

edge_range

2-tuple – Energy range (start, stop) used to determine the official beginning and edge of the edge itself.

E_0 = None

all_energies ()

edge_range = None

```
energies_in_range (norm_range=None)
map_range = None
post_edge = None
post_edge_order = 2
pre_edge = None
pre_edge_fit = None
regions = []
```

class `xanespy.edges.KEdge`
Bases: `xanespy.edges.Edge`
An X-ray absorption K-edge corresponding to a 1s transition.

annotate_spectrum (*ax*)

mask (**args, **kwargs*)
Return a numpy array mask for material that's active at this edge. Calculations are done in `xanes_math.l_edge_mask()`.

class `xanespy.edges.LEdge`
Bases: `xanespy.edges.Edge`
An X-ray absorption K-edge corresponding to a 2s or 2p transition.

annotate_spectrum (*ax*)

mask (**args, **kwargs*)
Return a numpy array mask for material that's active at this edge. Calculations are done in `xanes_math.l_edge_mask()`.

class `xanespy.edges.LMOMnKEdge`
Bases: `xanespy.edges.KEdge`
regions = [(6450, 6510, 20), (6524, 6542, 2), (6544, 6564, 1), (6566, 6568, 2), (6572, 6600, 4), (6610, 6650, 10), (6700, 6850,

class `xanespy.edges.NCACobaltLEdge`
Bases: `xanespy.edges.LEdge`
E_0 = 793.2
map_range = (0, 1)
post_edge = (785, 790)
pre_edge = (770, 775)
regions = [(770, 775, 1), (775, 785, 0.5), (785, 790, 1)]

class `xanespy.edges.NCANickelKEdge`
Bases: `xanespy.edges.KEdge`
E_0 = 8333
edge_range = (8341, 8358)
map_range = (8341, 8358)
post_edge = (8440, 8640)
pre_edge = (8249, 8281)
regions = [(8250, 8310, 20), (8324, 8344, 2), (8344, 8356, 1), (8356, 8360, 2), (8360, 8400, 4), (8400, 8440, 8), (8440, 8640, 5,

```

    shell = 'K'

class xanespy.edges.NCANickelKEdge61
    Bases: xanespy.edges.NCANickelKEdge

    regions = [(8250, 8310, 15), (8324, 8360, 1), (8360, 8400, 4), (8400, 8440, 8), (8440, 8640, 50)]

class xanespy.edges.NCANickelKEdge62
    Bases: xanespy.edges.NCANickelKEdge

    regions = [(8250, 8310, 15), (8324, 8360, 1), (8360, 8400, 4), (8400, 8440, 8), (8440, 8690, 50)]

class xanespy.edges.NCANickelLEdge
    Bases: xanespy.edges.LEdge

    E_0 = 853

    map_range = (0, 1)

    post_edge = (857, 862)

    pre_edge = (844, 848)

    regions = [(844, 848, 1), (849, 856, 0.25), (857, 862, 1)]

```

6.1.4 xanespy.exceptions module

Define classes for more fine-grained control over exception handling.

exception `xanespy.exceptions.CreateGroupError`

Bases: `ValueError`

Tried to import a TXM frameset into a group but the corresponding HDF group already exists or is otherwise inaccessible.

exception `xanespy.exceptions.DataFormatError`

Bases: `RuntimeError`

The raw data are arranged in a way that the importers or TXM classes do not understand.

exception `xanespy.exceptions.DataNotFoundError`

Bases: `FileNotFoundError`

Expected a directory containing data but found none.

exception `xanespy.exceptions.DatasetExistsError`

Bases: `RuntimeError`

Trying to save a new dataset but one already exists with the given path.

exception `xanespy.exceptions.FileExistsError`

Bases: `OSError`

Tried to import a TXM frameset but the corresponding HDF file already exists.

exception `xanespy.exceptions.FilenameParseError`

Bases: `ValueError`

The parameters in the filename do not match the naming scheme associated with this flavor.

exception `xanespy.exceptions.FrameFileNotFound`

Bases: `OSError`

Expected to load a TXM frame file but it doesn't exist.

exception `xanespy.exceptions.GroupKeyError`

Bases: `KeyError`

Tried to load or create an HDF group but failed. Examples include: the group doesn't exist, is ambiguous or already exists when being created.

exception `xanespy.exceptions.HDFScopeError`

Bases: `ValueError`

Tried to pass an HDF scope that is not recognized.

exception `xanespy.exceptions.NoParticleError`

Bases: `Exception`

exception `xanespy.exceptions.RefinementError`

Bases: `RuntimeError`

exception `xanespy.exceptions.XanesMathError`

Bases: `RuntimeError`

6.1.5 xanespy.importers module

`xanespy.importers.decode_aps_params(filename)`

Accept the filename of an XRM file and return sample parameters as a dictionary.

`xanespy.importers.decode_ssrl_params(filename)`

Accept the filename of an XRM file and return sample parameters as a dictionary.

`xanespy.importers.import_aps_8BM_frameset(directory, hdf_filename, quiet=False)`

`xanespy.importers.import_frameset(directory, flavor, hdf_filename)`

Import all files in the given directory collected at an X-ray microscope beamline.

Images are assumed to full-field transmission X-ray micrographs.

Parameters

- **directory** (*str*) – A valid path to a directory containing the frame data to import.
- **flavor** (*str*) – Indicates what type of naming conventions and data structure to assume. See documentation for `xanespy.xradia.XRMFile` for possible choice.
- **hdf_filename** (*str*) – Where to save the output to. An exception is throw if this file already exists.

`xanespy.importers.import_nanosurveyor_frameset(directory: str, quiet=False, hdf_filename=None, hdf_groupname=None, energy_range=None, exclude_re=None, append=False)`

Import a set of images from reconstructed ptychography scanning microscope data.

This generates ptychography chemical maps based on data collected at ALS beamline 5.3.2.1

Parameters

- **directory** (*str*) – Directory where to look for results. It should contain .cxi files that are the output of the ptychography reconstruction.”
- **quiet** (*Bool, optional*) – If truthy, progress bars will not be shown.
- **hdf_filename** (*str, optional*) – HDF File used to store computed results. If omitted or None, the *directory* basename is used

- **hdf_groupname** (*str*, *optional*) – Name to use for the hdf group of this dataset. If omitted or None, the *directory* basename is used. Raises an exception if the group already exists in the HDF file.
- **energy_range** (*2-tuple*, *optional*) – A 2-tuple with the (min, max) energy to be imported. This is useful if only a subset of the available data is usable. Values are assumed to be in electron-volts.
- **exclude_re** (*str*, *optional*) – Any filenames matching this regular expression will not be imported. A string or compiled re object can be given.
- **append** (*bool*, *optional*) – If True, any existing dataset will be added to, rather than replaced (default False)

`xanespy.importers.import_ssrl_frameset` (*directory*, *hdf_filename=None*)

Import all files in the given directory collected at SSRL beamline 6-2c and process into framesets. Images are assumed to full-field transmission X-ray micrographs and repetitions will be averaged. Passed on to `xanespy.importers.import_frameset`

`xanespy.importers.magnification_correction` (*frames*, *pixel_sizes*)

Correct for changes in magnification at different energies.

As the X-ray energy increases, the focal length of the zone plate changes and so the image is zoomed-out at higher energies. This method applies a correction to each frame to make the magnification similar to that of the first frame. Some beamlines correct for this automatically during acquisition and don't need this function: APS 8-BM-B, 32-ID-C.

Parameters

- **frames** (*np.ndarray*) – Numpy array of image frames that need to be corrected.
- **pixel_sizes** (*np.ndarray*) – Numpy array of pixel sizes corresponding to entries in *frames*.

Returns (*scales2D*, *translations*) – An array of scale factors to use for applying a correction to each frame. Translations show how much to move each frame array by to re-center it.

Return type (*np.ndarray*, *np.ndarray*)

`xanespy.importers.read_metadata` (*filenames*, *flavor*)

Take a list of filenames and return a pandas dataframe with all the metadata.

6.1.6 xanespy.plots module

Helper functions for setting up and displaying plots using matplotlib.

`xanespy.plots.big_axes` ()

Return a new Axes object, but larger than the default.

`xanespy.plots.draw_colorbar` (*ax*, *cmap*, *norm*, *energies*, *orientation='vertical'*, **args*, ***kwargs*)

Draw a colorbar on the side of a mapping axes to show the range of colors used. Returns the newly created colorbar object.

Parameters

- **ax** (–) –
- **cmap** (–) – to use.
- **norm** (–) – use.
- **energies** (–) – colorbar.

`xanespy.plots.draw_histogram_colorbar` (*ax*, *args, **kwargs)

Similar to `draw_colorbar()` with some special formatting options to put it along the X-axis of the axes.

`xanespy.plots.dual_axes` (*fig=None*, *orientation='horizontal'*)

Two new axes for mapping, side-by-side.

`xanespy.plots.new_axes` (*height=5*, *width=None*)

Create a new set of matplotlib axes for plotting. Height in inches.

`xanespy.plots.new_image_axes` (*height=5*, *width=5*)

Square axes with ticks on the outside.

`xanespy.plots.plot_composite_map` (*data*, *ax=None*, *origin='lower'*, *args, **kwargs)

Plot an RGB composite map on the given axes.

`xanespy.plots.plot_pixel_spectra` (*pixels*, *extent*, *spectra*, *energies*, *map_ax*, *spectra_ax=None*,
step_size=0)

Highlight certain pixels in an already-plotted map and plot their spectra. The map should already have been plotted.

Parameters

- **pixels** (-) – to highlight and plot.
- **extent** (-) – positions to (row, column) positions.
- **spectra** (-) – in *pixels* will use this array to get spectra. Shape is assumed to be (row, column, energy).
- **energies** (-) – plotting spectra.
- **map_ax** (-) – pixels.
- **spectra_ax** (-) – None (default) a new axes will be created.
- **step_size** (-) – directly on top of each other.

`xanespy.plots.plot_txm_histogram` (*data*, *ax=None*, *norm=None*, *bins=None*, *cmap='plasma'*,
add_cbar=True, *args, **kwargs)

Take an array of data values and show a histogram with some color-coding related to normalization value.

Returns: The matplotlib axes object used for plotting.

Parameters

- **data** (*np.ndarray*) – An array of values to plot on the histogram.
- **ax** (*optional*) – Matplotlib Axes instance to receive the plot. If None, a new axes will be created.
- **norm** (*optional*) – Matplotlib Normalize instance with the colormap range.
- **bins** (*optional*) – Bins to pass to the matplotlib `hist()` routine. If None (default), we will choose based on dtype of the data: integers will yield 1-wide bins, anything else will give 256 bins.
- **cmap** (*str, optional*) – Matplotlib colormap for coloring the bars.
- **add_cbar** (*bool, optional*) – Boolean to decide whether to add a colorbar along the bottom axis or not.
- ***args** – Positional arguments passed to matplotlib's `hist` call.
- ***kwargs** – Keyword arguments passed to matplotlib's `hist` call.

`xanespy.plots.plot_txm_intermediates` (*images*)

Accept a dictionary of images and plots them each on its own axes using matplotlib's *imshow*. This is a complement to routines that operate on a microscopy frame and optionally return all the intermediate calculated frames.

`xanespy.plots.plot_txm_map` (*data*, *edge=None*, *norm=None*, *ax=None*, *cmap='plasma'*, *origin='upper'*, **args*, ***kwargs*)

`xanespy.plots.plot_xanes_spectrum` (*spectrum*, *energies*, *norm=<matplotlib.colors.Normalize object>*, *show_fit=False*, *ax=None*, *ax2=None*, *linestyle=':'*, *color='blue'*, *cmap='plasma'*)

Plot a XANES spectrum on an axes. Applies some color formatting if *edge* is a valid XANES Edge object.

Parameters

- **spectrum** (–) –
- **energies** (–) –
- **norm** (–) – range. This will be used to annotate the plot if it is give.
- **show_fit** (–) –
- **ax** (–) – will be generated.
- **ax2** (–) – the data are complex.
- **linestyle** (–) –
- **cmap** (–) –
- **color** (–) – or “y” will decide based on the numerical value, *norm* and *cmap* arguments. Anything else will be passed as a color spec to the matplotlib commands.

`xanespy.plots.remove_extra_spines` (*ax*)

Removes the right and top borders from the axes.

`xanespy.plots.set_axes_color` (*ax*, *color*)

Set the axes, tick marks, etc of *ax* to mpl color *color*. Also, “doegreen” has special significance as the color associated with the US department of energy.

`xanespy.plots.set_outside_ticks` (*ax*)

Convert all the axes so that the ticks are on the outside and don't obscure data.

6.1.7 xanespy.qt_frame_view module

class `xanespy.qt_frame_view.FrameAnimation` (*fig*, *artists*, **args*, ***kwargs*)

Bases: `matplotlib.animation.ArtistAnimation`

Performs the animation for scrolling through frames arbitrarily.

stop ()

class `xanespy.qt_frame_view.FrameChangeSource` (*view*, **args*, ***kwargs*)

Bases: `PyQt5.QtCore.QObject`

add_callback (*func*, **args*, ***kwargs*)

callbacks = []

remove_callback (*func*, **args*, ***kwargs*)

start ()

stop ()

```
class xanespy.qt_frame_view.QtFrameView
    Bases: PyQt5.QtCore.QObject

    add_hdf_tree_item(item)

    clear_axes()

    connect_signals(presenter)

    create_canvas()

    create_status_bar()

    disable_frame_controls(status)

    disable_plotting_controls(status)

    draw_frames

    draw_histogram

    draw_spectrum(spectrum, energies, norm, cmap, edge_range)

    edge_ax = None

    expand_hdf_tree

    fig = None

    frame_changed

    frame_controls
        Gives a list of all the UI buttons that are associated with changing the currently active frame.

    hist_ax = None

    hist_cb = None

    plotting_controls
        Gives a list of all the UI elements that are associated with changing how the frameset is plotted.

    select_active_hdf_item(item)

    set_cmap(cmap)

    set_cmap_list(cmap_list)

    set_drawing_status(status)

    set_slider_max(val)

    set_status_cursor(msg)

    set_status_energy(msg)

    set_status_index(msg)

    set_status_pixel(msg)

    set_status_shape(msg)

    set_status_unit(msg)

    set_status_value(msg)

    set_timestep(idx)

    set_timestep_list(timestep_list)
```

set_ui_enabled (*enable=True*)
Turn on (default) or off the main interactive elements in the frame window. Useful for indicating blocking operations.

set_vmax (*val*)

set_vmax_decimals (*val*)

set_vmax_minimum (*val*)

set_vmax_step (*val*)

set_vmin (*val*)

set_vmin_decimals (*val*)

set_vmin_maximum (*val*)

set_vmin_step (*val*)

set_window_title (*title*)

setup ()

show ()

show_status_message (*message*)

spectrum_ax = None

ui = None

use_busy_cursor (*status*)

window = None

6.1.8 xanespy.qt_frameset_presenter module

class xanespy.qt_frameset_presenter.**QtFramesetPresenter** (*frameset, frame_view, *args, **kwargs*)

Bases: PyQt5.QtCore.QObject

Presenter for showing XanesFrameset frames and maps via Qt.

app_ready

pyqtSignal – Emitted when the application has been created and is ready for drawing.

busy_status_changed

pyqtSignal – Emitted when processing has started or ended.

map_data_changed

pyqtSignal – Emitted when the map data is different and should be re-plotted.

map_data_cleared

pyqtSignal – Emitted when no map data is available and plots can be cleared.

active_frame = 0

active_frames ()

active_map ()

Returns the active map array if possible. If it doesn't exist, return None.

active_representation = None

active_timestep = 0

add_frame_view (*view*, *threaded=True*)

Attach a view to this presenter.

Parameters

- **view** (*QObject*) – The view will be connected to signals that describe changes in frame data.
- **threaded** (*bool*, *optional*) – If true, this view will be added to its own thread before signals get connected

add_map_view (*view*, *threaded=True*)

Attach a view to this presenter.

Parameters

- **view** (*QObject*) – The view will be connected to signals that describe changes in map data.
- **threaded** (*bool*, *optional*) – If true, this view will be added to its own thread before signals get connected, giving a snappier UI. Disabling makes testing more straightforward.

animate_frames ()

app_ready

build_hdf_tree ()

Build the items and insert them into the view’s HDF tree based on the structure of the frameset’s HDF file.

busy_status_changed

change_cmap (*new_cmap*)

change_hdf_group (*new_item*, *old_item*)

change_map_cmap (*new_cmap*)

cmap_list_changed

connect_signals ()

create_app ()

draw_frame_histogram ()

draw_frame_spectra ()

first_frame ()

frame_cmap = ‘copper’

frame_norm ()

frame_pixel = None

hover_frame_pixel (*xy*)

last_frame ()

launch ()

map_cmap = ‘plasma’

map_cursor_changed

map_data_changed

map_data_cleared

map_limits_changed

map_norm()

map_pixel_changed

map_spectrum_changed

mean_spectrum_changed

move_map_pixel (*vert*, *horiz*)

Move the active pixel by the given amount in each direction.

If the current pixel is not active, this is a no-op.

move_slider (*new_value*)

next_frame()

num_frames = 0

play_frames (*start*)

prepare_ui()

previous_frame()

process_events

refresh_frames()

reset_frame_range()

Reset the frame plotting vmin and vmax based on the currently selected data. VMin will be the 1st percentile and VMax will be the 99th percentile. The results can be accessed by calling the *frame_norm()* method.

Returns

- (**vmin**, **vmax**) ((*float*, *float*)) – A 2-tuple with the new min and max range.
- *None* – If current representation is not a valid map, None will be returned.

reset_map_range()

Reset the map plotting vmin and vmax based on the currently selected data. VMin will be the 1st percentile and VMax will be the 99th percentile. This method also caches the values so they can be retrieved via the *map_norm()* method. If the current representation is not valid map data, then the cached vmin and vmax will be returned with no other changes.

Returns (**vmin**, **vmax**) – A 2-tuple with the current min and max range.

Return type (*float*, *float*)

set_frame_vmax (*new_value*)

set_frame_vmin (*new_value*)

set_map_cursor (*x*, *y*)

set_map_pixel (*x*, *y*)

set_map_vmax (*new_value*)

set_map_vmin (*new_value*)

set_play_speed (*new_speed*)

Change how fast the play timer ticks. Input speeds should be in the range of 0, 30 with 30 being the fastest. This is converted to timer intervals between 1ms and 1000ms on an exponential scale.

set_timestep (*new_timestep*)

timer = None

toggle_edge_mask (*state*)

update_frame_range_limits ()

Check that the (min, max, step) of the frame spinboxes are reasonable. This should be called when the spinbox values change.

update_map_limits ()

update_maps ()

Send the current mapping data to the *map_data_changed* signal.

This method should be called after anything changes the visual representation of the map. Examples:

- Changing the active representation
- Changing the colormap
- Changing the map normalization limits
- Changing whether the XAS edge mask is applied

update_spectra ()

Get the most recent data for mean and single-pixel spectra and send them out to the signals *mean_spectrum_changed* and *pixel_spectrum_changed*.

update_status_frame (*new_frame*)

Create a string (and send it to the UI) that indicates the energy of the requested frame.

update_status_shape ()

update_status_unit ()

update_status_value ()

use_edge_mask = False

6.1.9 xanespy.qt_map_view module

class xanespy.qt_map_view.**QtMapView**

Bases: PyQt5.QtCore.QObject

A Qt view for a frameset map. It should be controlled by a presenter.

cmap_changed

signal – Fires when the user changes the colormap via the UI

limits_applied

signal – Fires when the user requests that data be redrawn with new limits.

limits_reset

signal – Fires when the user asks that the norm limits be reset to the data.

cmap_changed

connect_presenter (*presenter*)

Connect to signals for changed presenter state.

create_canvas ()

crosshairs = None

edge_mask_toggled

```

fig = None
hide()
keyboard_nav(event)
latest_cmap = 'plasma'
limits_applied
limits_reset
map_clicked
map_hovered
map_moved
map_vmax_changed
map_vmin_changed
mouse_clicked_canvas(mouse_event)
mouse_in_canvas(mouse_event)
plot_histogram_data(map_data, norm, cmap, extent)
plot_map_data(map_data, norm, cmap, extent)
plot_spectrum(spectrum, norm, cmap, edge_range)
redraw_canvas()
redraw_crosshairs(xy)
    Draw a set of crosshairs on the map at location given by xy.
set_cmap_list(new_list)
set_map_limits(vmin, vmax, step, decimals)
setup_ui()
show()
ui = None
update_crosshair_labels(xy, pixel, value)
update_cursor_labels(xy, pixel, value)
window = None

```

6.1.10 xanespy.txmstore module

Tools for accessing TXM data stored in an HDF5 file.

```

class xanespy.txmstore.TXMStore(hdf_filename: str, parent_name: str, data_name=None,
                                mode='r')
    Bases: object

```

Wrapper around HDF5 file that stores TXM data. It has a series of properties that return the corresponding HDF5 dataset object; the `TXMStore().attribute.value` pattern can be used to get pure numpy arrays. These objects should be used as a context manager to ensure that the file is closed, especially if using a writing mode:

```

with TXMStore() as store: # Do stuff with store here

```

Parameters

- **hdf_filename** (*str*) – Path to the HDF file to be used.
- **parent_name** (*str*) – Name of the top-level HDF5 group.
- **data_name** (*str*) – Name of the second level HDF5 group, used for specific data iterations (eg. imported, aligned)
- **mode** (*str*) – Eg. ‘r’ for read-only, ‘r+’ for read-write. Passed directly to h5py.File constructor.

VERSION = 1

absorbances

close()

cluster_map

data_group()

Retrieve the currently active second-level HDF5 group object for this file and groupname. Ex. “imported” or “aligned_frames”.

data_name

data_tree()

Create a tree of the possible groups this store could access. The first level is samples, then data_groups (ie. same sample but different analysis status), then representations. Maps are not included in this tree.

energies

filenames

fit_parameters

fork_data_group (*new_name*)

Turn on different active data group for this store. This method deletes the existing group and copies symlinks from the current one.

get_dataset (*name*)

Attempt to open the requested dataset.

Returns **data** – An open HDF5 dataset

Return type h5py.Dataset

Raises `exceptions.GroupKeyError` – If the dataset does not exist in the file.

get_frames (*name*)

Get a set of frames, specified by the value of *name*.

get_map (*name*)

Get a map of the frames, specified by the value of *name*.

has_dataset (*name*)

Return a boolean indicated whether this dataset exists in the HDF file.

intensities

latest_data_name

original_positions

parent_group()

Retrieve the top-level HDF5 group object for this file and groupname.

particle_labels**pixel_sizes****pixel_unit****references****relative_positions**

(x, y, z) position values for each frame.

replace_dataset (*name, data, context=None, attrs={}, compression=None, *args, **kwargs*)Wrapper for `h5py.create_dataset` that removes the existing dataset if it exists.**Parameters**

- **name** (*str*) – HDF5 groupname name to give this dataset.
- **data** (*np.ndarray*) – Numpy array of data to be saved.
- **context** (*str, optional*) – Specifies what kind of data is stored. Eg. “frameset”, “metadata”, “map”.
- **attrs** (*dict, optional*) – Dictionary containing HDF5 metadata attributes to be set on the resulting dataset.
- ***args** – Arguments to pass to `h5py’s create_dataset` method.
- ****kwargs** – Keyword arguments to pass to `h5py’s create_dataset` method.

set_frames (*name, val*)Set data for a set of frames, specified by the value of *name*.**signal_map****signal_method**

String describing how the previously extracted signals were calculated.

signal_weightsGet the pixel weights of the previously extracted signals using any one of a variety of decomposition methods, saved as *signal_method*.**signals**Get the previously extracted signals using any one of a variety of decomposition methods, saved as *signal_method*.**timestamps****timestep_names****whiteline_map**

6.1.11 xanespy.utilities module

A collection of classes and functions that aren’t specific to any one type of measurement. Also, it defines some namedtuples for describing coordinates.

xanespy.utilities.Extent

alias of extent

class xanespy.utilities.Pixel (*vertical, horizontal*)

Bases: tuple

horizontal

Alias for field number 1

vertical

Alias for field number 0

`xanespy.utilities.broadcast_reverse` (*array, shape, *args, **kwargs*)

Take the array and extends it as much as possible to match *shape*. Similar to numpy's `broadcast_to` function, but starts with the most significant axis. For example, if *array* has shape (7, 29), it can be broadcast to (7, 29, 1024, 1024).

`xanespy.utilities.foreach` (*f, l, threads=6, return_=False*)

Apply *f* to each element of *l*, in parallel

`xanespy.utilities.get_component` (*data, name*)

If complex, turn to given component, otherwise return original data.

`xanespy.utilities.is_kernel` ()

Detect whether or not we're running inside an IPython kernel. NB: This does not distinguish between eg IPython notebook and IPython QtConsole.

`xanespy.utilities.parallel_map` (*f, l, threads=6*)

`xanespy.utilities.pixel_to_xy` (*pixel, extent, shape*)

Take an xy location on an image and convert it to a pixel location suitable for numpy indexing.

class `xanespy.utilities.position` (*x, y, z*)

Bases: tuple

x

Alias for field number 0

y

Alias for field number 1

z

Alias for field number 2

`xanespy.utilities.prog` (*leave=False, dynamic_ncols=True, *args, **kwargs*)

A progress bar for displaying how many iterations have been completed. This is mostly just a wrapper around the `tqdm` library. *args* and *kwargs* are passed directly to either `tqdm.tqdm` or `tqdm.tqdm_notebook`.

Parameters

- **leave** (*bool, optional*) – Whether to leave the progress bar in the stream after it's completed.
- **dynamic_ncols** (*bool, optional*) – Whether to adapt dynamically to the width of the environment.
- **args** – Positional arguments passed directly to `tqdm` or `tqdm_notebook`.
- **kwargs** – Keyword arguments passed directly to `tqdm` or `tqdm_notebook`.

class `xanespy.utilities.shape` (*rows, columns*)

Bases: tuple

columns

Alias for field number 1

rows

Alias for field number 0

`xanespy.utilities.xy_to_pixel` (*xy, extent, shape*)

Take an xy location on an image and convert it to a pixel location suitable for numpy indexing.

```
class xanespy.utilities.xycoord(x, y)
    Bases: tuple

    x
        Alias for field number 0

    y
        Alias for field number 1
```

6.1.12 xanespy.xanes_frameset module

Class definitions for working with a whole stack of X-ray microscopy frames. Each frame is a micrograph at a different energy. A frameset then is a three-dimensional dataset with of dimensions (energy, row, column).

```
class xanespy.xanes_frameset.PtychoFrameset(filename, edge, groupname=None)
    Bases: xanespy.xanes_frameset.XanesFrameset
```

A set of images (“frames”) at different energies moving across an absorption edge. The individual frames should be generated by ptychographic reconstruction of scanning transmission X-ray microscopy (STXM) to produce an array complex intensity values. This class does *not* include any code responsible for the collection and reconstruction of such data, only for the analysis in the context of X-ray absorption near edge spectroscopy.

```
apply_internal_reference(plot_background=True, ax=None)
```

Use a portion of each frame for internal reference correction. The result is the complex refraction for each pixel: the real component describes the phase shift, and the imaginary component is exponential decay, ie. absorbance.

Parameters

- **plot_background** (If *truthy*, the values of *I₀* are plotted as) – a function of energy.
- **ax** (The axes to use for plotting if *plot_background* is) – *truthy*.

```
representations()
```

Retrieve a list of valid representations for these data, such as modulus or phase data for ptychography.

```
class xanespy.xanes_frameset.XanesFrameset(filename, edge, groupname=None)
    Bases: object
```

A collection of TXM frames at different energies moving across an absorption edge. Iterating over this object gives the individual Frame() objects. The class assumes that the data have been imported into an HDF file.

```
active_group = “
```

```
active_labels_groupname
```

The group name for the latest frameset of detected particle labels.

```
add_frame(frame)
```

```
align_frames(reference_frame='mean', blur=None, method: str = 'cross_correlation', tem-
               plate=None, passes=1, commit=True, component='modulus', plot_results=True)
```

Use cross correlation algorithm to line up the frames. All frames will have their sample position set to (0, 0) since we don’t know which one is the real position. This operation will interpolate between pixels so introduces error. If multiple passes are performed, the translations are saved and combined at the end so this error is only introduced once. Using the *commit=False* argument allows for multiple different types of registration to be performed in sequence, since uncommitted translations will be applied before the next round of registration.

Parameters

- **(int, str or None)** (*reference_frame*) – which all other frames should be aligned. If None, the frame of highest intensity will be used. If “mean” (default) or “median”, the average or median of all frames will be used. If “max”, the frame with highest absorbance is used. This attribute has no effect if template matching is used.
- **blur** (*A type of filter to apply to each frame of the data*) – before attempting registration. Choices are “median” or None
- **method** (*Which technique to use to calculate the translation*) –
 - “cross_correlation” (default)
 - “template_match”(If “template_match” is used, the *template* argument should also be provided.)
- **passes** (*How many times this alignment should be done. Default: 1.*) –
- **template** (*Image data that should be matched if the*) – *template_match* method is used.
- **commit** (*If truthy (default), the final translation will be*) – applied to the data stored on disk by calling *self.apply_translations(crop=True)* after all passes have finished.
- **component** (*What component of the data to use: 'modulus',*) – ‘phase’, ‘imag’ or ‘real’.
- **plot_results** (*If truthy (default), plot the root-mean-square of the*) – translation distance for each pass.

apply_transformations (*crop=True, commit=True*)

Take any transformations staged with *self.stage_transformations()* and apply them. If commit is truthy, the staged transformations are reset.

Returns: Transformed array of the absorbances frames.

Parameters

- **crop** (–) –
- **so there are not edges. If falsy, the images will** (*translated,*) –
- **wrapped.** (*be*) –
- **commit** (–) – store for absorbances, intensities and references, and the staged transformations will be cleared. Otherwise, only the absorbance data will be transformed and returned.
- **frames_name** (–) – transformation too (eg. ‘absorbances’)

background_group ()

background_normalizer ()

calculate_maps ()

Generate a set of maps based on pixel-wise Xanes spectra: whiteline position, particle labels.

Parameters **-fit_spectra** (*If truthy, the whiteline will be found by*)
– fitting curves, instead of the default of taking the direct maximum.

calculate_signals (*n_components=2, method='nmf', edge_mask=True*)

Extract signals and assign each pixel to a group, then save the resulting RGB cluster map.

Parameters

- **n_components** (-) – into which the data will be separated.
- **method** (-) – signals. Currently only “nmf” is supported.
- **edge_mask** (-) – the edge filter will be considered.

calculate_whitelines (*edge_mask=False*)

Calculate and save a map of the whiteline position of each pixel by calculating the energy of simple maximum absorbance.

Parameters **edge_mask** (-) – be fit and the remaning pixels will be set to a default value. This can help reduce computing time.

clear_caches ()

Clear cached function values so they will be recomputed with fresh data

cmap = ‘plasma’

data_name

data_tree ()

Wrapper around the TXMStore.data_tree() method.

drop_frame (*index*)

Delete frame with the given index (int) or energy (float). This destructively removes the data from the HDF5 file, so use with caution.

edge_jump

Calculate a image showing the difference in signal across the X-ray edge.

edge_mask

Calculate a mask for what is likely active material at this edge.

Parameters

- **sensitivity** (*float*) – A multiplier for the otsu value to determine the actual threshold.
- **min_size** (*int*) – Objects below this size (in pixels) will be removed. Passing zero (default) will result in no effect.

endtime ()

Determine the latest timestamp amongst all of the frames.

energies

Return the array of beam energies for the given time index.

Returns **energies** – A 1-dimensional array with the energy for each frame.

Return type np.ndarray

extent

Determine physical dimensions for axes values.

Returns: If idx was given, a single tuple of (left, right, bottom, up), otherwise if idx is None it returns an array of extents for each frame.

Parameters **-idx** (*Index for a given frame. This allows for faster*) – calculation if only a single frame is required. By default, returns the first frame.

fit_spectra (*edge_mask=True*)

Fit a series of curves to the spectrum at each pixel.

Parameters `edge_mask` (*bool, optional*) – If true, only pixels passing the `edge_mask` will be fit and the remaining pixels will be set to a default value. This can help reduce computing time.

fork_data_group (*new_name, src=None*)

Turn on different active data for this frameset's store object. Similar to `switch_data_group` except that this method deletes the existing group and copies symlinks from the current one.

Parameters

- **new_name** (-) –
- **src** (-) – from. If None (default), the current data group will be used.

frame_shape (*representation='absorbances'*)

Return the shape of the individual energy frames.

frames

Return the frames for the given time index.

If *representation* is really mapping data, then the source frames will be returned.

Parameters

- **timeidx** (*int*) – Index for the first dimension of the combined data array.
- **representation** (*str*) – The group name for these data. Eg “absorbances”, “white-line_map”, “intensities”

Returns **frames** – A 3-dimensional array with the form (energy, row, col).

Return type np.ndarray

goodness_mask (*sensitivity: float = 0.5*)

Calculate an image based on the goodness of fit. `calculate_map` must have been called previously. Goodness filter is converted to a binary map using (Otsu) threshold filtering.

Parameters **sensitivity** (-) – the actual threshold.

gtk_viewer ()

Launch a GTK gui to view the data in the frameset.

has_representation (*representation*)

hdf_file (*mode='r'*)

Return an open h5py.File object for this (any maybe other) frameset.

Parameters

- **mode** (A mode string, see h5py documentation for options. To)–
- **file corruption, calling this method as a context** (*avoid*)–
- **is recommended, especially with a mode other than 'r'**.
(*manager*) –

hdf_node ()

For use with HDFAttribute descriptor.

hdf_path (*representation=None*)

Return the hdf path for the active group.

Parameters **representation** (*str, optional*) – Name of third-level group to use. If omitted, the path to the parent group will be given.

Returns `path` – The path to the current group in the HDF5 file. Returns an empty string if the representation does not exists.

Return type `str`

image_normalizer

is_background()

label_particles (*min_distance=20*)

Use watershed segmentation to identify particles.

Parameters `min_distance` (–) – grouping areas into particles. Lower numbers means more particles, but might split large particles into two.

map_data

Return map data for the given time index and representation.

If *representation* is really mapping data, then the result will have more dimensions than expected.

Parameters

- **timeidx** (*int*) – Index for the first dimension of the combined data array.
- **representation** (*str*) – The group name for these data. Eg “absorbances”, “white-line_map”, “intensities”

Returns `map_data` – A 2-dimensional array with the form (row, col).

Return type `np.ndarray`

masked_map (*goodness_filter=True*)

Generate a map based on pixel-wise Xanes spectra and apply an edge-jump filter mask. Default is to compute X-ray whiteness position.

mean_frame()

Return the mean absorbance with the same shape as an individual frame.

movie_plotter()

Creates an animation of all the frames in ascending energy, but does not display it anywhere, that’s up to you.

normalize (*plot_fit=False, new_name='normalized'*)

Correct for background material not absorbing at this edge. Uses method described in DOI 10.1038/ncomms7883: fit line against material that fails `edge_jump_filter` and use this line to correct entire frame.

Parameters

- **plot_fit** (–) –
- **line.** (*best-fit*) –

num_energies

num_timesteps

particle (*particle_idx=0*)

Prepare a particle frameset for the given particle index.

particle_area_spectrum (*loc=xycoord(x=20, y=20)*)

Calculate a spectrum based on the area of the particle closest to the given location in the frame. This may be useful for assessing magnification across multiple frames.

particle_centroid_spectrum (*loc=xycoord(x=20, y=20)*)

Calculate a spectrum based on the image centroid of the particle closest to the given location in the frame. This may be useful for assessing systematic drift across multiple frames.

particle_regions (*intensity_image=None, labels=None*)

Return a list of regions (1 for each particle) sorted by area. (largest first). This requires that the *label_particles* method be called first.

Parameters

- **intensity_image** (-) – *regionprops* function to determine what shows up in the image for each particle.
- **labels** (-) – particles segmented. If None, the *particle_labels* attribute of the TXM store will be used.

particle_series (*map_name='whiteline_map'*)

Generate an array of values from *map_name* averaged across each particle.

Returns: A 2D array where the first dimension is particles and the second is the first dimension of the map dataset (usually time).

pixel_unit ()

Return the unit of measure for the size of a pixel.

plot_edge_jump (*ax=None, alpha=1*)

Plot the results of the edge jump filter.

plot_frame (*idx, ax=None, cmap='gray', *args, **kwargs*)

Plot the frame with given index as an image.

plot_goodness (*plotter=None, ax=None, norm_range=None, *args, **kwargs*)

Use a default frameset plotter to draw a map of the goodness of fit as determined by the Edge object.

plot_histogram (*plotter=None, timeidx=None, ax=None, vmin=None, vmax=None, goodness_filter=False, active_pixel=None, bins='energies', *args, **kwargs*)

Use a default frameset plotter to draw a map of the chemical data.

plot_map (*ax=None, map_name='whiteline_map', timeidx=0, vmin=None, vmax=None*)

Prepare data and plot a map of whiteline positions.

plot_map_pixel_spectra (*pixels, map_ax=None, spectra_ax=None, map_name='whiteline_map', timeidx=0, step_size=0*)

Plot the frameset's map and highlight some pixels on it then plot those pixel's spectra on another set of axes.

Parameters

- **pixels** (-) – column pixels to highlight.
- **map_ax** (-) – None, a new 2-wide subplot will be created for both *map_ax* and *spectra_ax*.
- **spectra_ax** (-) – spectra. Will only be used if *map_ax* is not None.
- **map_name** (-) – passed to the TXM store object and retrieved from the hdf5 file. If falsy, no map will be plotted.
- **timeidx** (-) –

plot_mean_image (*ax=None, component='modulus'*)

plot_signal_map (*ax=None, signals_idx=None*)

Plot the map of signal strength for signals extracted from *self.calculate_signals()*.

Parameters

- **ax** (-) – axes object is created.
- **signals_idx** (-) –
- **as a numpy array index. Special value None(default) (passed)** –
- **first three signals will be plotted. (means)** –

plot_signals (*cmap='viridis'*)

Plot the signals from the previously extracted data. Requires that `self.store().signals` and `self.store().signal_weights` be set.

plot_xanes_edge (**args, **kwargs*)

Call `self.plots.plot_xanes_spectrum()` but zoomed in on the edge.

plot_xanes_spectrum (*ax=None, pixel=None, norm_range=None, normalize=False, representation='modulus', show_fit=False, edge_jump_filter=False, linestyle=':', *args, **kwargs*)

Calculate and plot the xanes spectrum for this field-of-view.

Parameters

- **- matplotlib axes object on which to draw (ax)** –
- **- Coordinates of a specific pixel on the image to plot. (pixel)** –
- **- If truthy, will set the pre-edge at zero and the (normalize)** – post-edge at 1.
- **- If truthy, will use the edge object to fit the data (show_fit)** – and plot the resulting fit line.
- **- If truthy, will only include those values (edge_jump_filter)** – that show a strong absorbtion jump across this edge.

qt_viewer ()

representations ()

Retrieve a list of valid representations for these data.

save_movie (*filename, *args, **kwargs*)

Save an animation of all the frames and XANES to the specified filename.

spectra (*edge_filter=False*)

Return a two-dimensional array of spectra for all the pixels in shape of (pixel, energy).

Parameters (bool or str) (edge_jump_filter) – truthy, only pixels that pass the edge jump filter are used to calculate the spectrum. If “inverse” is given, then the edge jump filter is logically not-ted and calculated with a more conservative threshold.

spectrum (*pixel=None, edge_jump_filter=False, representation='absorbances', index=0*)

Collapse the frameset down to an energy spectrum.

Any dimensions (besides the energy dimension) that remain after applying the arguments below, will be averaged to give the final intensity at each energy.

Returns spectrum – A pandas Series with the spectrum.

Return type `pd.Series`

Parameters

- **pixel** (*tuple, optional*) – A 2-tuple that causes the returned series to represent the spectrum for only 1 pixel in the frameset. If None, a larger part of the frame will be used, depending on the other arguments.
- **edge_jump_filter** (*bool or str, optional*) – If truthy, only pixels that pass the edge jump filter are used to calculate the spectrum. If “inverse” is given, then the edge jump filter is logically not-ted and calculated with a more conservative threshold.
- **representation** (*str, optional*) – What kind of data to use for creating the spectrum. This will be passed to `TXMStore.get_map()`
- **index** (*int, optional*) – Which step in the frameset to use. When used to index `store().absorbances`, this should return a 3D array like (energy, rows, columns).

stage_transformations (*translations=None, rotations=None, center=(0, 0), scales=None*)

Allows for deferred transformation of the frame data.

Since each transformation introduces interpolation error, the best results occur when the translations are saved up and then applied all in one shot. Takes a combination of arrays of translations (x, y), rotations and/or scales and saves them for later application. This method should be used in conjunction `apply_transformations()`.

All three arguments should have shapes that are compatible with the frame data, though this is not strictly enforced for now. Rotation will necessarily have one less degree of freedom than translation/scale values.

Example Shapes:

Frames	Translations	Rotations	Scales
(10, 48, 1024, 1024)	(10, 48, 2)	(10, 48, 1)	(10, 48, 2)
(10, 48, 1024, 1024, 1024)	(10, 48, 3)	(10, 48, 2)	(10, 48, 3)

Parameters

- **translations** (*np.ndarray*) – How much to move each axis (x, y[, z]).
- **rotations** (*np.ndarray*) – How much to rotate around the origin (0, 0) pixel.
- **center** (*2-tuple*) – Where to set the origin of rotation. Default is the first pixel (0, 0).
- **scales** (*np.ndarray*) – How much to scale the image by in each dimension (x, y[, z]).

starttime ()

Determine the earliest timestamp amongst all of the frames.

store (*mode='r'*)

Get a TXM Store object that saves and retrieves data from the HDF5 file. The mode argument is passed to `h5py` as is. This method should be used as a context manager, especially if mode is something writeable:

with self.store() as store: # Do stuff with the store

subtract_surroundings ()

Use the edge mask to separate “surroundings” from “sample”, then subtract the average surrounding absorbance from each frame. This effectively removes effects where the entire frame is brighter from one energy to the next.

6.1.13 xanespy.xanes_math module

Module containing all the computationally demanding functions. This allows for easy optimization of parallelizable algorithms. Most functions will operate on large arrays of data.

class xanespy.xanes_math.**KEdgeParams** (*scale, voffset, E0, sigw, pre_m, pre_b, ga, gb, gc*)

Bases: tuple

E0

Alias for field number 2

ga

Alias for field number 6

gb

Alias for field number 7

gc

Alias for field number 8

pre_b

Alias for field number 5

pre_m

Alias for field number 4

scale

Alias for field number 0

sigw

Alias for field number 3

voffset

Alias for field number 1

xanespy.xanes_math.**apply_internal_reference** (*intensities, out=None*)

Apply a reference correction to complex data to convert intensities into refractive index. I_0 is determined by separating the pixels into background and foreground using Otsu's method.

Arrays *intensities* and *out* must all have the same shape where the last two dimensions are image rows and column.

xanespy.xanes_math.**apply_references** (*intensities, references, out=None*)

Apply a reference correction to convert intensity values to optical depth.

The formula $-\ln \frac{\text{intensities}}{\text{references}}$ is used to calculate the new values. Arrays *intensities*, *references* and *out* must all have the same shape where the last two dimensions are image rows and columns.

Parameters

- **intensities** (*np.ndarray*) – Sample input signal data.
- **references** (*np.ndarray*) – Background input signal data. Must be the same shape as *intensities*.
- **out** (*np.ndarray, optional*) – Array to receive the results.

xanespy.xanes_math.**direct_whitelines** (*spectra, energies, edge*)

Takes an array of X-ray absorbance spectra and calculates the positions of maximum intensities over the near-edge region.

Parameters

- **spectra** (*np.array*) – 2D numpy array of absorbance spectra where the last dimension is energy.
- **energies** (*np.array*) – Array of X-ray energies in electron-volts. Must be broadcastable to the shape of *spectra*.
- **edge** – An XAS Edge object that describes the absorbance edge in question.

Returns **out** – Array with the whitenline position of each spectrum.

Return type np.ndarray

`xanespy.xanes_math.extract_signals_nmf(spectra, n_components, nmf_kwargs=None, mask=None)`

Extract the signal components present in the given spectra using non-negative matrix factorization. Input data can be negative, but it will be shifted up, processed, then shifted down again.

Parameters

- **spectra** (–) –
- **n_components** (–) –
- **nmf_kwargs** (–) – the constructor of the estimator.

Returns

Return type 2-tuple of arrays (components, weights)

`xanespy.xanes_math.fit_kedge(spectra, energies, p0)`

Use least squares to fit a set of curves to the data. Currently this is a line for the baseline absorbance decreasing at higher energies, plus a sigmoid for the edge and a gaussian for the whitenline.

Returns an array with a similar shape to spectra but the last axis is replaced with fitting parameters, describe by the named tuple *KParams* defined in this module.

Parameters

- **spectra** (–) – index is energy. This can be a multi-dimensional array, which allows calculation of image frames, etc. The last axis should be X-ray energy.
- **energies** (–) –
- **p0** (–) – described by *kedge_params*.
- **out** (–) –
- **be created.** (*will*) –

`xanespy.xanes_math.guess_kedge(spectrum, energies, edge)`

Guess initial starting parameters for a k-edge curve. This will give a rough estimate, appropriate for giving to the *fit_kedge* function as the starting parameters, *p0*.

Parameters

- **spectrum** (–) – K-edge spectrum. Only 1-dimensional data are currently accepted.
- **energies** (–) – the points in *spectrum*. Must have the same shape as *spectrum*.
- **edge** (–) – actual edge energy itself.
- **Returns** (A named tuple with the estimated parameters (see) – *.KEdgeParams* for definition)

`xanespy.xanes_math.iter_indices(data, leftover_dims=1, desc=None)`

Accept an array of frames, indices, etc. and generate slices for each frame. Assumes the last two dimensions of *data* are rows and columns. All other dimensions will be iterated over.

- **leftover_dims** : Integer describing which dimensions should not be iterated over. Eg. if data is 3D array and *leftover_dims* == 1, only first two dimensions will be iterated.
- **desc** : String to put in the progress bar.

`xanespy.xanes_math.k_edge_jump(frames: numpy.ndarray, energies: numpy.ndarray, edge)`

Determine what the difference is between the *post_edge* and the *pre_edge*.

`xanespy.xanes_math.k_edge_mask` (*frames: numpy.ndarray, energies: numpy.ndarray, edge, sensitivity: float = 1, min_size=0*)

Calculate a mask for what is likely active material at this edge. This is done by comparing the edge-jump to the standard deviation. Foreground material will be identified when the edge-jump accounts for most of the standard deviation.

Parameters

- **frames** (-) –
- **energies** (-) – *frames*. Must have the same shape along the first dimension as *frames*.
- **edge** (-) – elemental edge being studied.
- **sensitivity** (-) – the actual threshold.
- **min_size** (-) – removed. Passing zero (default) will result in no effect.

Returns

- - A boolean mask with the same shape as the last two dimensions of
- *frames* where True pixels are likely to be background material.

`xanespy.xanes_math.l_edge_mask` (*frames: numpy.ndarray, energies: numpy.ndarray, edge, sensitivity: float = 1, frame_dims=2, min_size=0*)

Calculate a mask for what is likely active material at this edge. This is done by comparing each spectrum to the overall spectrum using the dot product. A normalization is first applied to mitigate differences in total intensity.

Parameters

- **frames** (-) –
- **energies** (-) – *frames*. Must have the same shape along the first dimension as *frames*.
- **edge** (-) – elemental edge being studied.
- **sensitivity** (-) – the actual threshold.
- **frame_dims** (-) – means each frame is a two-dimensional image.
- **min_size** (-) – removed. Passing zero (default) will result in no effect.

Returns

- - A boolean mask with the same shape as the last two dimensions of
- *frames* where True pixels are likely to be background material.

`xanespy.xanes_math.particle_labels` (*frames: numpy.ndarray, energies: numpy.ndarray, edge, min_distance=20*)

Prepare a map by segmenting the images into particles.

Parameters **frames** (-) – energy. These will be merged and used for segmentation.

`xanespy.xanes_math.predict_edge` (*energies, *params*)

Defines the curve function that gets fit to the data for an absorbance K-edge.

The predicted curve is a combination of a straight line (for background), an arctan function (for the edge), and a gaussian peak (for the whitenline).

Parameters

- **energies** (*np.ndarray*) – Array with energy values to be predicted.
- ***params** (*tuple(int)*) – The curve parameters that should be used for the prediction. Their order is described by `kedge_params` variable.

Returns *curve* – The predicted absorbance values based on the input parameters. Shape will match *energies*.

Return type `np.ndarray`

`xanespy.xanes_math.register_correlations` (*frames*, *reference*, *upsample_factor*=10, *desc*='Registering')

Calculate the relative translation between the reference image and a series of frames.

This uses phase correlation through scikit-image's *register_translation* function.

Parameters

- **frames** (*np.ndarray*) – Array where the last two dimensions are (column, row) of images to be registered.
- **reference** (*np.ndarray*) – Image frame against which to align the entries in *frames*.
- **upsample_factor** (*int*, *optional*) – Factor controls subpixel registration via scikit-image.
- **desc** (*str*, *optional*) – Description for putting in the progress bar.

Returns *translations* – Array with same dimensions as 0-th axis of *frames* containing (x, y) translations for each frame.

Return type `np.ndarray`

`xanespy.xanes_math.register_template` (*frames*, *reference*, *template*, *desc*='Registering')

Calculate the relative translation between the reference image and a series of frames.

This uses template cross correlation through scikit-image's *match_template* function.

The *register_correlations* algorithm is simpler to use in most cases but sometimes results in unreasonable results; in those cases, this method can be more reliable to achieve a first approximation.

Parameters

- **frames** (*np.ndarray*) – Array where the last two dimensions are (column, row) of images to be registered.
- **reference** (*np.ndarray*) – Image frame against which to align the entries in *frames*.
- **template** (*np.ndarray*) – A 2D array (smaller than frames and reference) that will be identified in each frame and used for alignment.
- **desc** (*str*, *optional*) – Description for putting in the progress bar.

Returns *translations* – Array with same dimensions as 0-th axis of *frames* containing (x, y) translations for each frame.

Return type `np.ndarray`

`xanespy.xanes_math.transform_images` (*data*, *transformations*, *out*=None, *mode*='median')

Takes image data and applies the given translation matrices.

It is assumed that the first dimension of *data* is the same as the length of *transformations*. The transformation matrices can be generated from translation, rotation and scale parameters via the `xanespy.xanes_math.transformation_matrices()` function. Data will be written to *out* if given, otherwise returned as a new array.

Parameters

- **data** (*np.ndarray*) – Numeric array with frames to transform. Last two dimensions are assumed to be (row, columns).

- **transformations** (*np.ndarray*) – A numeric array shaped compatibly with *data*. The last two dimensions are assumed to be (3, 3) and each (3, 3) encodes a transformation matrix for the corresponding frame in *data*.
- **out** (*np.ndarray*, *optional*) – A numeric array with same shape as *data* that will hold the transformed data.
- **mode** (*str*, *optional*) – Describes how to deal with edges. See scikit-image documentation for options. Special value “median” (default), takes the median pixel intensity of that frame and uses it as the constant value.

Returns out – A new array with similar dimensions to *data* but with transformations applied and converted to float datatype.

Return type *np.ndarray*

`xanespy.xanes_math.transformation_matrices` (*translations=None*, *rotations=None*, *scales=None*, *center=(0, 0)*)

Takes array of operations and calculates (3, 3) transformation matrices.

This function operates by calculating an AffineTransform similar to that described in the scikit-image package.

All three arguments (*translations*, *rotations*, and *scales*, should have shapes that are compatible with the frame data, though this is not strictly enforced for now. Rotation will necessarily have one less degree of freedom than translation/scale values.

Example Shapes:

Frames	Translations	Rotations	Scales
(10, 48, 1024, 1024)	(10, 48, 2)	(10, 48, 1)	(10, 48, 2)
(10, 48, 1024, 1024, 1024)	(10, 48, 3)	(10, 48, 2)	(10, 48, 3)

Parameters

- **translations** (*np.ndarray*, *optional*) – How much to move each axis (x, y[, z]).
- **rotations** (*np.ndarray*, *optional*) – How much to rotate around the origin (0, 0) pixel.
- **center** (*np.ndarray*, *optional*) – Where to set the origin of rotation. Default is the first pixel (0, 0).
- **scales** (*np.ndarray*, *optional*) – How much to scale the image by in each dimension (x, y[, z]).

Returns new_transforms – Resulting transformation matrices. Will have the same shape as the input arrays but with the last dimension replaced by (3, 3).

Return type *np.ndarray*

6.1.14 xanespy.xradia module

Tools for importing X-ray microscopy frames in formats produced by Xradia instruments.

class `xanespy.xradia.XRMFile` (*filename*, *flavor: str*)

Bases: `object`

Single X-ray microscopy frame created using XRadia XRM format.

Parameters

- **filename** (–) –

- **flavor** (-) – choices are ['ssrl', 'aps', 'aps-old1']. These choices should line up with whatever is generated using the scripts in beamlines moudles.

aps_old1_regex = `re.compile('(\d{8})_([a-zA-Z0-9]+)_([a-zA-Z0-9]+)(\d{4}).xrm')`

binning ()

close ()

Close original XRM (ole) file on disk.

endtime ()

Retrieve a datetime object representing when this frame was finished collecting. Duration is decided by exposure time of the frame and the start time.

energy ()

Beam energy in electronvoltes.

image_data ()

TXM Image frame.

image_dtype ()

image_shape ()

ole_value (*stream, fmt=None*)

Get arbitrary data from the ole file and convert from bytes.

print_ole ()

sample_position ()

starttime ()

Retrieve a datetime object representing when this frame was collected. Timezone is inferred from flavor (eg. ssrl -> california time).

um_per_pixel ()

Describe the size of a pixel in microns. If this is an SSRL frame, the pixel size is dependent on energy. For APS frames, the pixel size is uniform and assumes a 40µm field-of-view.

6.1.15 Module contents

INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

PYTHON MODULE INDEX

X

- `xanespy`, [44](#)
- `xanespy.beamlines`, [13](#)
- `xanespy.edges`, [15](#)
- `xanespy.exceptions`, [17](#)
- `xanespy.importers`, [18](#)
- `xanespy.plots`, [19](#)
- `xanespy.qt_frame_view`, [21](#)
- `xanespy.qt_frameset_presenter`, [23](#)
- `xanespy.qt_map_view`, [26](#)
- `xanespy.txmstore`, [27](#)
- `xanespy.utilities`, [29](#)
- `xanespy.xanes_frameset`, [31](#)
- `xanespy.xanes_math`, [38](#)
- `xanespy.xradia`, [43](#)

INDEX

A

absorbances (xanespy.txmstore.TXMStore attribute), 28

active_frame (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 23

active_frames() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 23

active_group (xanespy.xanes_frameset.XanesFrameset attribute), 31

active_labels_groupname (xanespy.xanes_frameset.XanesFrameset attribute), 31

active_map() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 23

active_representation (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 23

active_timestep (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 23

add_callback() (xanespy.qt_frame_view.FrameChangeSource method), 21

add_frame() (xanespy.xanes_frameset.XanesFrameset method), 31

add_frame_view() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 23

add_hdf_tree_item() (xanespy.qt_frame_view.QtFrameView method), 22

add_map_view() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24

align_frames() (xanespy.xanes_frameset.XanesFrameset method), 31

all_energies() (xanespy.edges.Edge method), 15

animate_frames() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24

annotate_spectrum() (xanespy.edges.KEdge method), 16

annotate_spectrum() (xanespy.edges.LEdge method), 16

app_ready (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 23, 24

apply_internal_reference() (in module xane-

spy.xanes_math), 39

apply_internal_reference() (xanespy.xanes_frameset.PtychoFrameset method), 31

apply_references() (in module xanespy.xanes_math), 39

apply_transformations() (xanespy.xanes_frameset.XanesFrameset method), 32

aps_old1_regex (xanespy.xradia.XRMFile attribute), 44

B

background_group() (xanespy.xanes_frameset.XanesFrameset method), 32

background_normalizer() (xanespy.xanes_frameset.XanesFrameset method), 32

big_axes() (in module xanespy.plots), 19

binning() (xanespy.xradia.XRMFile method), 44

broadcast_reverse() (in module xanespy.utilities), 30

build_hdf_tree() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24

busy_status_changed (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 23, 24

C

calculate_maps() (xanespy.xanes_frameset.XanesFrameset method), 32

calculate_signals() (xanespy.xanes_frameset.XanesFrameset method), 32

calculate_whitelines() (xanespy.xanes_frameset.XanesFrameset method), 33

callbacks (xanespy.qt_frame_view.FrameChangeSource attribute), 21

change_cmap() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24

- change_hdf_group() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24
 change_map_cmap() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24
 clear_axes() (xanespy.qt_frame_view.QtFrameView method), 22
 clear_caches() (xanespy.xanes_frameset.XanesFrameset method), 33
 close() (xanespy.txmstore.TXMStore method), 28
 close() (xanespy.xradia.XRMFile method), 44
 cluster_map (xanespy.txmstore.TXMStore attribute), 28
 cmap (xanespy.xanes_frameset.XanesFrameset attribute), 33
 cmap_changed (xanespy.qt_map_view.QtMapView attribute), 26
 cmap_list_changed (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 24
 columns (xanespy.utilities.shape attribute), 30
 connect_presenter() (xanespy.qt_map_view.QtMapView method), 26
 connect_signals() (xanespy.qt_frame_view.QtFrameView method), 22
 connect_signals() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24
 create_app() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24
 create_canvas() (xanespy.qt_frame_view.QtFrameView method), 22
 create_canvas() (xanespy.qt_map_view.QtMapView method), 26
 create_status_bar() (xanespy.qt_frame_view.QtFrameView method), 22
 CreateGroupError, 17
 crosshairs (xanespy.qt_map_view.QtMapView attribute), 26
- ## D
- data_group() (xanespy.txmstore.TXMStore method), 28
 data_name (xanespy.txmstore.TXMStore attribute), 28
 data_name (xanespy.xanes_frameset.XanesFrameset attribute), 33
 data_tree() (xanespy.txmstore.TXMStore method), 28
 data_tree() (xanespy.xanes_frameset.XanesFrameset method), 33
 DataFormatError, 17
 DataNotFoundError, 17
 DatasetExistsError, 17
 decode_aps_params() (in module xanespy.importers), 18
 decode_ssrl_params() (in module xanespy.importers), 18
 Detector (class in xanespy.beamlines), 13
 DetectorPoint (class in xanespy.beamlines), 13
 direct_whitelines() (in module xanespy.xanes_math), 39
 disable_frame_controls() (xanespy.qt_frame_view.QtFrameView method), 22
 disable_plotting_controls() (xanespy.qt_frame_view.QtFrameView method), 22
 draw_colorbar() (in module xanespy.plots), 19
 draw_frame_histogram() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24
 draw_frame_spectra() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24
 draw_frames (xanespy.qt_frame_view.QtFrameView attribute), 22
 draw_histogram (xanespy.qt_frame_view.QtFrameView attribute), 22
 draw_histogram_colorbar() (in module xanespy.plots), 19
 draw_spectrum() (xanespy.qt_frame_view.QtFrameView method), 22
 drop_frame() (xanespy.xanes_frameset.XanesFrameset method), 33
 dual_axes() (in module xanespy.plots), 20
- ## E
- E_0 (xanespy.xanes_math.KEdgeParams attribute), 39
 E_0 (xanespy.edges.Edge attribute), 15
 E_0 (xanespy.edges.NCACobaltLEdge attribute), 16
 E_0 (xanespy.edges.NCANickelKEdge attribute), 16
 E_0 (xanespy.edges.NCANickelLEdge attribute), 17
 Edge (class in xanespy.edges), 15
 edge_ax (xanespy.qt_frame_view.QtFrameView attribute), 22
 edge_jump (xanespy.xanes_frameset.XanesFrameset attribute), 33
 edge_mask (xanespy.xanes_frameset.XanesFrameset attribute), 33
 edge_mask_toggled (xanespy.qt_map_view.QtMapView attribute), 26
 edge_range (xanespy.edges.Edge attribute), 15
 edge_range (xanespy.edges.NCANickelKEdge attribute), 16
 endtime() (xanespy.xanes_frameset.XanesFrameset method), 33
 endtime() (xanespy.xradia.XRMFile method), 44
 energies (xanespy.txmstore.TXMStore attribute), 28
 energies (xanespy.xanes_frameset.XanesFrameset attribute), 33
 energies_in_range() (xanespy.edges.Edge method), 16
 energy (xanespy.beamlines.DetectorPoint attribute), 13
 energy (xanespy.beamlines.ZoneplatePoint attribute), 13

energy() (xanespy.xradia.XRMFile method), 44
expand_hdf_tree (xanespy.qt_frame_view.QtFrameView attribute), 22
Extent (in module xanespy.utilities), 29
extent (xanespy.xanes_frameset.XanesFrameset attribute), 33
extract_signals_nmf() (in module xanespy.xanes_math), 40

F

fig (xanespy.qt_frame_view.QtFrameView attribute), 22
fig (xanespy.qt_map_view.QtMapView attribute), 27
FileExistsError, 17
FilenameParseError, 17
filenames (xanespy.txmstore.TXMStore attribute), 28
first_frame() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24
fit_kedge() (in module xanespy.xanes_math), 40
fit_parameters (xanespy.txmstore.TXMStore attribute), 28
fit_spectra() (xanespy.xanes_frameset.XanesFrameset method), 33
foreach() (in module xanespy.utilities), 30
fork_data_group() (xanespy.txmstore.TXMStore method), 28
fork_data_group() (xanespy.xanes_frameset.XanesFrameset method), 34
frame_changed (xanespy.qt_frame_view.QtFrameView attribute), 22
frame_cmap (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 24
frame_controls (xanespy.qt_frame_view.QtFrameView attribute), 22
frame_norm() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24
frame_pixel (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 24
frame_shape() (xanespy.xanes_frameset.XanesFrameset method), 34
FrameAnimation (class in xanespy.qt_frame_view), 21
FrameChangeSource (class in xanespy.qt_frame_view), 21
FrameFileNotFound, 17
frames (xanespy.xanes_frameset.XanesFrameset attribute), 34

G

ga (xanespy.xanes_math.KEdgeParams attribute), 39
gb (xanespy.xanes_math.KEdgeParams attribute), 39
gc (xanespy.xanes_math.KEdgeParams attribute), 39
get_component() (in module xanespy.utilities), 30
get_dataset() (xanespy.txmstore.TXMStore method), 28
get_frames() (xanespy.txmstore.TXMStore method), 28

get_map() (xanespy.txmstore.TXMStore method), 28
goodness_mask() (xanespy.xanes_frameset.XanesFrameset method), 34
GroupKeyError, 17
gtk_viewer() (xanespy.xanes_frameset.XanesFrameset method), 34
guess_kedge() (in module xanespy.xanes_math), 40

H

has_dataset() (xanespy.txmstore.TXMStore method), 28
has_representation() (xanespy.xanes_frameset.XanesFrameset method), 34
hdf_file() (xanespy.xanes_frameset.XanesFrameset method), 34
hdf_node() (xanespy.xanes_frameset.XanesFrameset method), 34
hdf_path() (xanespy.xanes_frameset.XanesFrameset method), 34
HDFScopeError, 18
hide() (xanespy.qt_map_view.QtMapView method), 27
hist_ax (xanespy.qt_frame_view.QtFrameView attribute), 22
hist_cb (xanespy.qt_frame_view.QtFrameView attribute), 22
horizontal (xanespy.utilities.Pixel attribute), 29
hover_frame_pixel() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24

I

image_data() (xanespy.xradia.XRMFile method), 44
image_dtype() (xanespy.xradia.XRMFile method), 44
image_normalizer (xanespy.xanes_frameset.XanesFrameset attribute), 35
image_shape() (xanespy.xradia.XRMFile method), 44
import_aps_8BM_frameset() (in module xanespy.importers), 18
import_frameset() (in module xanespy.importers), 18
import_nanosurveyor_frameset() (in module xanespy.importers), 18
import_sslr_frameset() (in module xanespy.importers), 19
intensities (xanespy.txmstore.TXMStore attribute), 28
is_background() (xanespy.xanes_frameset.XanesFrameset method), 35
is_kernel() (in module xanespy.utilities), 30
iter_indices() (in module xanespy.xanes_math), 40

K

k_edge_jump() (in module xanespy.xanes_math), 40
k_edge_mask() (in module xanespy.xanes_math), 40

KEdge (class in xanespy.edges), 16

KEdgeParams (class in xanespy.xanes_math), 38

keyboard_nav() (xanespy.qt_map_view.QtMapView method), 27

L

l_edge_mask() (in module xanespy.xanes_math), 41

label_particles() (xanespy.xanes_frameset.XanesFrameset method), 35

last_frame() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24

latest_cmap (xanespy.qt_map_view.QtMapView attribute), 27

latest_data_name (xanespy.txmstore.TXMStore attribute), 28

launch() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 24

LEdge (class in xanespy.edges), 16

limits_applied (xanespy.qt_map_view.QtMapView attribute), 26, 27

limits_reset (xanespy.qt_map_view.QtMapView attribute), 26, 27

LMOMnKEdge (class in xanespy.edges), 16

M

magnification_correction() (in module xanespy.importers), 19

map_clicked (xanespy.qt_map_view.QtMapView attribute), 27

map_cmap (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 24

map_cursor_changed (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 24

map_data (xanespy.xanes_frameset.XanesFrameset attribute), 35

map_data_changed (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 23, 24

map_data_cleared (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 23, 24

map_hovered (xanespy.qt_map_view.QtMapView attribute), 27

map_limits_changed (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 24

map_moved (xanespy.qt_map_view.QtMapView attribute), 27

map_norm() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25

map_pixel_changed (xanespy.qt_frameset_presenter.QtFramesetPresenter

attribute), 25

map_range (xanespy.edges.Edge attribute), 15, 16

map_range (xanespy.edges.NCACobaltLEdge attribute), 16

map_range (xanespy.edges.NCANickelKEdge attribute), 16

map_range (xanespy.edges.NCANickelLEdge attribute), 17

map_spectrum_changed (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 25

map_vmax_changed (xanespy.qt_map_view.QtMapView attribute), 27

map_vmin_changed (xanespy.qt_map_view.QtMapView attribute), 27

mask() (xanespy.edges.KEdge method), 16

mask() (xanespy.edges.LEdge method), 16

masked_map() (xanespy.xanes_frameset.XanesFrameset method), 35

mean_frame() (xanespy.xanes_frameset.XanesFrameset method), 35

mean_spectrum_changed (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 25

mouse_clicked_canvas() (xanespy.qt_map_view.QtMapView method), 27

mouse_in_canvas() (xanespy.qt_map_view.QtMapView method), 27

move_map_pixel() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25

move_slider() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25

movie_plotter() (xanespy.xanes_frameset.XanesFrameset method), 35

N

name (xanespy.edges.Edge attribute), 15

NCACobaltLEdge (class in xanespy.edges), 16

NCANickelKEdge (class in xanespy.edges), 16

NCANickelKEdge61 (class in xanespy.edges), 17

NCANickelKEdge62 (class in xanespy.edges), 17

NCANickelLEdge (class in xanespy.edges), 17

new_axes() (in module xanespy.plots), 20

new_image_axes() (in module xanespy.plots), 20

next_frame() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25

NoParticleError, 18

normalize() (xanespy.xanes_frameset.XanesFrameset method), 35

num_energies (xanespy.xanes_frameset.XanesFrameset attribute), 35

num_frames (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 25
 num_timesteps (xanespy.xanes_frameset.XanesFrameset attribute), 35

O

ole_value() (xanespy.xradia.XRMFile method), 44
 original_positions (xanespy.txmstore.TXMStore attribute), 28

P

parallel_map() (in module xanespy.utilities), 30
 parent_group() (xanespy.txmstore.TXMStore method), 28
 particle() (xanespy.xanes_frameset.XanesFrameset method), 35
 particle_area_spectrum() (xanespy.xanes_frameset.XanesFrameset method), 35
 particle_centroid_spectrum() (xanespy.xanes_frameset.XanesFrameset method), 35
 particle_labels (xanespy.txmstore.TXMStore attribute), 28
 particle_labels() (in module xanespy.xanes_math), 41
 particle_regions() (xanespy.xanes_frameset.XanesFrameset method), 36
 particle_series() (xanespy.xanes_frameset.XanesFrameset method), 36
 Pixel (class in xanespy.utilities), 29
 pixel_sizes (xanespy.txmstore.TXMStore attribute), 29
 pixel_to_xy() (in module xanespy.utilities), 30
 pixel_unit (xanespy.txmstore.TXMStore attribute), 29
 pixel_unit() (xanespy.xanes_frameset.XanesFrameset method), 36
 play_frames() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25
 plot_composite_map() (in module xanespy.plots), 20
 plot_edge_jump() (xanespy.xanes_frameset.XanesFrameset method), 36
 plot_frame() (xanespy.xanes_frameset.XanesFrameset method), 36
 plot_goodness() (xanespy.xanes_frameset.XanesFrameset method), 36
 plot_histogram() (xanespy.xanes_frameset.XanesFrameset method), 36
 plot_histogram_data() (xanespy.qt_map_view.QtMapView method), 27
 plot_map() (xanespy.xanes_frameset.XanesFrameset method), 36
 plot_map_data() (xanespy.qt_map_view.QtMapView method), 27
 plot_map_pixel_spectra() (xanespy.xanes_frameset.XanesFrameset method), 36
 plot_mean_image() (xanespy.xanes_frameset.XanesFrameset method), 36
 plot_pixel_spectra() (in module xanespy.plots), 20
 plot_signal_map() (xanespy.xanes_frameset.XanesFrameset method), 36
 plot_signals() (xanespy.xanes_frameset.XanesFrameset method), 37
 plot_spectrum() (xanespy.qt_map_view.QtMapView method), 27
 plot_txm_histogram() (in module xanespy.plots), 20
 plot_txm_intermediates() (in module xanespy.plots), 20
 plot_txm_map() (in module xanespy.plots), 21
 plot_xanes_edge() (xanespy.xanes_frameset.XanesFrameset method), 37
 plot_xanes_spectrum() (in module xanespy.plots), 21
 plot_xanes_spectrum() (xanespy.xanes_frameset.XanesFrameset method), 37
 plotting_controls (xanespy.qt_frame_view.QtFrameView attribute), 22
 position (class in xanespy.utilities), 30
 position() (xanespy.beamlines.Zoneplate method), 13
 post_edge (xanespy.edges.Edge attribute), 15, 16
 post_edge (xanespy.edges.NCACobaltLEdge attribute), 16
 post_edge (xanespy.edges.NCANickelKEdge attribute), 16
 post_edge (xanespy.edges.NCANickelLEdge attribute), 17
 post_edge_order (xanespy.edges.Edge attribute), 15, 16
 pre_b (xanespy.xanes_math.KEdgeParams attribute), 39
 pre_edge (xanespy.edges.Edge attribute), 15, 16
 pre_edge (xanespy.edges.NCACobaltLEdge attribute), 16
 pre_edge (xanespy.edges.NCANickelKEdge attribute), 16
 pre_edge (xanespy.edges.NCANickelLEdge attribute), 17
 pre_edge_fit (xanespy.edges.Edge attribute), 16
 pre_m (xanespy.xanes_math.KEdgeParams attribute), 39
 predict_edge() (in module xanespy.xanes_math), 41
 prepare_ui() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25
 previous_frame() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25
 print_ole() (xanespy.xradia.XRMFile method), 44
 process_events (xanespy.qt_frameset_presenter.QtFramesetPresenter

attribute), 25
 prog() (in module xanespy.utilities), 30
 PtychoFrameset (class in xanespy.xanes_frameset), 31

Q

qt_viewer() (xanespy.xanes_frameset.XanesFrameset method), 37
 QtFramesetPresenter (class in xanespy.qt_frameset_presenter), 23
 QtFrameView (class in xanespy.qt_frame_view), 21
 QtMapView (class in xanespy.qt_map_view), 26

R

read_metadata() (in module xanespy.importers), 19
 redraw_canvas() (xanespy.qt_map_view.QtMapView method), 27
 redraw_crosshairs() (xanespy.qt_map_view.QtMapView method), 27
 references (xanespy.txmstore.TXMStore attribute), 29
 RefinementError, 18
 refresh_frames() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25
 regions (xanespy.edges.Edge attribute), 15, 16
 regions (xanespy.edges.LMOMnKEdge attribute), 16
 regions (xanespy.edges.NCACobaltLEdge attribute), 16
 regions (xanespy.edges.NCANickelKEdge attribute), 16
 regions (xanespy.edges.NCANickelKEdge61 attribute), 17
 regions (xanespy.edges.NCANickelKEdge62 attribute), 17
 regions (xanespy.edges.NCANickelLEdge attribute), 17
 register_correlations() (in module xanespy.xanes_math), 42
 register_template() (in module xanespy.xanes_math), 42
 relative_positions (xanespy.txmstore.TXMStore attribute), 29
 remove_callback() (xanespy.qt_frame_view.FrameChangeSource method), 21
 remove_extra_spines() (in module xanespy.plots), 21
 replace_dataset() (xanespy.txmstore.TXMStore method), 29
 representations() (xanespy.xanes_frameset.PtychoFrameset method), 31
 representations() (xanespy.xanes_frameset.XanesFrameset method), 37
 reset_frame_range() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25
 reset_map_range() (xanespy.qt_frameset_presenter.QtFramesetPresenter

method), 25
 rows (xanespy.utilities.shape attribute), 30

S

sample_position() (xanespy.xradia.XRMFile method), 44
 save_movie() (xanespy.xanes_frameset.XanesFrameset method), 37
 scale (xanespy.xanes_math.KEdgeParams attribute), 39
 sector8_xanes_script() (in module xanespy.beamlines), 14
 select_active_hdf_item() (xanespy.qt_frame_view.QtFrameView method), 22
 set_axes_color() (in module xanespy.plots), 21
 set_cmap() (xanespy.qt_frame_view.QtFrameView method), 22
 set_cmap_list() (xanespy.qt_frame_view.QtFrameView method), 22
 set_cmap_list() (xanespy.qt_map_view.QtMapView method), 27
 set_drawing_status() (xanespy.qt_frame_view.QtFrameView method), 22
 set_frame_vmax() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25
 set_frame_vmin() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25
 set_frames() (xanespy.txmstore.TXMStore method), 29
 set_map_cursor() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25
 set_map_limits() (xanespy.qt_map_view.QtMapView method), 27
 set_map_pixel() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25
 set_map_vmax() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25
 set_map_vmin() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25
 set_outside_ticks() (in module xanespy.plots), 21
 set_play_speed() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 25
 set_slider_max() (xanespy.qt_frame_view.QtFrameView method), 22
 set_status_cursor() (xanespy.qt_frame_view.QtFrameView method), 22
 set_status_energy() (xanespy.qt_frame_view.QtFrameView method),

[22](#)
 set_status_index() (xanespy.qt_frame_view.QtFrameView method), [22](#)
 set_status_pixel() (xanespy.qt_frame_view.QtFrameView method), [22](#)
 set_status_shape() (xanespy.qt_frame_view.QtFrameView method), [22](#)
 set_status_unit() (xanespy.qt_frame_view.QtFrameView method), [22](#)
 set_status_value() (xanespy.qt_frame_view.QtFrameView method), [22](#)
 set_timestep() (xanespy.qt_frame_view.QtFrameView method), [22](#)
 set_timestep() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), [25](#)
 set_timestep_list() (xanespy.qt_frame_view.QtFrameView method), [22](#)
 set_ui_enabled() (xanespy.qt_frame_view.QtFrameView method), [22](#)
 set_vmax() (xanespy.qt_frame_view.QtFrameView method), [23](#)
 set_vmax_decimals() (xanespy.qt_frame_view.QtFrameView method), [23](#)
 set_vmax_minimum() (xanespy.qt_frame_view.QtFrameView method), [23](#)
 set_vmax_step() (xanespy.qt_frame_view.QtFrameView method), [23](#)
 set_vmin() (xanespy.qt_frame_view.QtFrameView method), [23](#)
 set_vmin_decimals() (xanespy.qt_frame_view.QtFrameView method), [23](#)
 set_vmin_maximum() (xanespy.qt_frame_view.QtFrameView method), [23](#)
 set_vmin_step() (xanespy.qt_frame_view.QtFrameView method), [23](#)
 set_window_title() (xanespy.qt_frame_view.QtFrameView method), [23](#)
 setup() (xanespy.qt_frame_view.QtFrameView method), [23](#)
 setup_ui() (xanespy.qt_map_view.QtMapView method), [27](#)
 shape (class in xanespy.utilities), [30](#)
 shell (xanespy.edges.NCANickelKEdge attribute), [16](#)
 show() (xanespy.qt_frame_view.QtFrameView method), [23](#)
 show() (xanespy.qt_map_view.QtMapView method), [27](#)
 show_status_message() (xanespy.qt_frame_view.QtFrameView method), [23](#)
 signal_map (xanespy.txmstore.TXMStore attribute), [29](#)
 signal_method (xanespy.txmstore.TXMStore attribute), [29](#)
 signal_weights (xanespy.txmstore.TXMStore attribute), [29](#)
 signals (xanespy.txmstore.TXMStore attribute), [29](#)
 sigw (xanespy.xanes_math.KEdgeParams attribute), [39](#)
 spectra() (xanespy.xanes_frameset.XanesFrameset method), [37](#)
 spectrum() (xanespy.xanes_frameset.XanesFrameset method), [37](#)
 spectrum_ax (xanespy.qt_frame_view.QtFrameView attribute), [23](#)
 ssl6_xanes_script() (in module xanespy.beamlines), [14](#)
 stage_transformations() (xanespy.xanes_frameset.XanesFrameset method), [38](#)
 start() (xanespy.qt_frame_view.FrameChangeSource method), [21](#)
 starttime() (xanespy.xanes_frameset.XanesFrameset method), [38](#)
 starttime() (xanespy.xradia.XRMFile method), [44](#)
 stop() (xanespy.qt_frame_view.FrameAnimation method), [21](#)
 stop() (xanespy.qt_frame_view.FrameChangeSource method), [21](#)
 store() (xanespy.xanes_frameset.XanesFrameset method), [38](#)
 subtract_surroundings() (xanespy.xanes_frameset.XanesFrameset method), [38](#)

T

timer (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), [26](#)
 timestamps (xanespy.txmstore.TXMStore attribute), [29](#)
 timestep_names (xanespy.txmstore.TXMStore attribute), [29](#)
 toggle_edge_mask() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), [26](#)
 transform_images() (in module xanespy.xanes_math), [42](#)
 transformation_matrices() (in module xanespy.xanes_math), [43](#)
 TXMStore (class in xanespy.txmstore), [27](#)

U

ui (xanespy.qt_frame_view.QtFrameView attribute), [23](#)
 ui (xanespy.qt_map_view.QtMapView attribute), [27](#)
 um_per_pixel() (xanespy.xradia.XRMFile method), [44](#)

update_crosshair_labels() (xanespy.qt_map_view.QtMapView method), 27
 update_cursor_labels() (xanespy.qt_map_view.QtMapView method), 27
 update_frame_range_limits() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 26
 update_map_limits() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 26
 update_maps() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 26
 update_spectra() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 26
 update_status_frame() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 26
 update_status_shape() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 26
 update_status_unit() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 26
 update_status_value() (xanespy.qt_frameset_presenter.QtFramesetPresenter method), 26
 use_busy_cursor() (xanespy.qt_frame_view.QtFrameView method), 23
 use_edge_mask (xanespy.qt_frameset_presenter.QtFramesetPresenter attribute), 26

V

VERSION (xanespy.txmstore.TXMStore attribute), 28
 vertical (xanespy.utilities.Pixel attribute), 29
 voffset (xanespy.xanes_math.KEdgeParams attribute), 39

W

whiteline_map (xanespy.txmstore.TXMStore attribute), 29
 window (xanespy.qt_frame_view.QtFrameView attribute), 23
 window (xanespy.qt_map_view.QtMapView attribute), 27
 write_scaninfo_header() (in module xanespy.beamlines), 15

X

x (xanespy.beamlines.ZoneplatePoint attribute), 14
 x (xanespy.utilities.position attribute), 30
 x (xanespy.utilities.xycoord attribute), 31
 XanesFrameset (class in xanespy.xanes_frameset), 31

XanesMathError, 18
 xanespy (module), 44
 xanespy.beamlines (module), 13
 xanespy.edges (module), 15
 xanespy.exceptions (module), 17
 xanespy.importers (module), 18
 xanespy.plots (module), 19
 xanespy.qt_frame_view (module), 21
 xanespy.qt_frameset_presenter (module), 23
 xanespy.qt_map_view (module), 26
 xanespy.txmstore (module), 27
 xanespy.utilities (module), 29
 xanespy.xanes_frameset (module), 31
 xanespy.xanes_math (module), 38
 xanespy.xradia (module), 43
 XRMFile (class in xanespy.xradia), 43
 xy_to_pixel() (in module xanespy.utilities), 30
 xycoord (class in xanespy.utilities), 30

Y

y (xanespy.beamlines.ZoneplatePoint attribute), 14
 y (xanespy.utilities.position attribute), 30
 y (xanespy.utilities.xycoord attribute), 31

Z

z (xanespy.beamlines.DetectorPoint attribute), 13
 z (xanespy.beamlines.ZoneplatePoint attribute), 14
 z (xanespy.utilities.position attribute), 30
 Zoneplate (class in xanespy.beamlines), 13
 ZoneplatePoint (class in xanespy.beamlines), 13